

BiPAP 101 For Non-respiratory Personnel

All the PAPs and what to do with them

Purpose

This training program is intended for non-respiratory personnel who work with patients who are on a Respirationics Vision BiPAP unit. It is intended to be very basic training and is in no way a training program that assures competence in using a BiPAP machine. Consult a trained respiratory therapist if there are any questions or problems.

The Three PAPs

BiPAP = Bi-level positive airway pressure or Non-invasive positive pressure ventilation (NPPV)

IPAP = Inspiratory Positive Airway Pressure. More IPAP (with same EPAP)= bigger breath = blow off more CO₂

EPAP = Expiratory Positive Airway Pressure. Same as PEEP or CPAP. More EPAP = better end-expiratory inflation = better PaO₂

Pressure support = IPAP – EPAP = breath boost = size of breath

The Masks

Total face mask

More comfortable than it looks

One size fits most

Avoid nasal pressure

Place on cannula at similar FiO₂ for eating etc

The Masks

Place strap behind head

Make sure unit is on or in the stand-by mode

The Masks

Place mask on face and support with one hand

The Masks

Adjust straps to seal leaks

Straps can be folded over to not interfere with vision

The Masks

Full face mask

Sized like an ambu bag mask

Watch for tissue break-down

Place on cannula at similar FiO₂ for eating etc.

The Masks

Nasal mask

Mostly for obstructive sleep apnea patients

Nasal pillows also used

Place on cannula at similar FiO₂ for eating etc

Mask Fitting (Parameters Screen)

0 – 6 liter per minute leak = Too tight (loosen mask)

6 – 30 lpm leak = just right

30 – 60 lpm leak = Acceptable but monitor volumes closely

> 60 lpm leak = Re-adjust mask

Putting on Stand-by

Used for short-term removals (eating, meds etc)

The Vision shows Stand-by on the right side of the MONITORING screen

Remove the mask from the patient

Press Standby

The machine will switch to BiPAP automatically when the mask is placed on the patient

Putting on Stand-by

Used for short-term removals (eating, meds etc)

The Vision shows Stand-by on the right side of the MONITORING screen

Remove the mask from the patient

Press Standby

The machine will switch to BiPAP automatically when the mask is placed on the patient

Putting on Stand-by

Place on cannula or mask at similar FiO₂ for eating etc.

Approximate FiO₂'s delivered by cannulas

1 liter = 24%

2 liters = 28%

3 liters = 32%

4 liters = 36%

5 liters = 40%

6 liters = 44%

Turning BiPAP Off

Disconnect patient

Push power button on back of machine

Call RT when machine needs to be turned on – they will do the pre-use check

Parameters

Changed and recorder by RT

Emergency F_iO₂ changes

Press parameters

Press FIO₂

Rotate knob

Notify RT immediately!

Managing Alarms

Audible alarm will sound and alarm message will be displayed at the top of the display screen

SILENCE mutes the alarm for 60 seconds

Reset clears alarm messages

Call RT to troubleshoot continuing alarms

Managing Alarms

Audible alarm will sound and alarm message will be displayed at the top of the display screen

SILENCE mutes the alarm for 60 seconds

Reset clears alarm messages

Call RT to troubleshoot continuing alarms

Managing Alarms – Vent Inop

Red wrench

Remove patient, assure patient safety and call RT

Managing Alarms – Check Ventilator

Yellow eye

Remove patient, assure patient safety and call RT

Conclusion

NPPV is safe, effective and widely used

Contact RT with any questions

Assure patient safety – BVM is as effective (short term)